

Genotyping with the Tag-It™ Mutation Detection Platform

Richard Janeczko, PhD
Salisbury, July 2004



MADRID TALK

- How to take penalty kicks

Overview

- **The Tag-It™ Genotyping* Platform**
- **Tag-It™ Product Menu**
 - Cystic fibrosis,
 - Ashkenazi Jewish
 - Thrombophilia
 - Drug metabolism (2D6, 2C9, 2C19)
- **In Development**
 - Fragile X
 - Infectious Disease (MRSA, viral genotypes)




Tag-It™ Product Specifications


- Accurate, robust, reproducible
- Flexible open system allows for
 - Easy addition and subtraction of specific mutations by testing lab
 - Open home-brew platform
- Variable sample throughput capability
- Software:
 - Reports both wild-type and mutant alleles
 - Integrates with LIMS and generates clinical reports
- Cost-effective




ASPE Genotyping with 'Tag-It' or how to use the Tag-It Array




Tag-It platform: product format



Tag-It, Universal Array Sorting



86 Luminex Bead Populations



TAG IT

Tag-It™ Cystic Fibrosis Products: Mutation Panels

- CFTR 23+4 (available March 5th 2004)
- CFTR 25+4
 - 25 most prevalent (ACMG set)+ 4 reflex tests, Q3 2003
- CFTR 40+4 set
 - ACMG set + 15 world's most common + 4 reflex tests, Q3 2003
- CFTR 71+4
 - CFTR 40+4 set + 31 expanded US and European coverage (in development), Q4 2004



Lab Requirements

4 Major Steps

Schematic

Hands On Time

I Multiplex PCR

Genomic DNA (25 ng)



SAP-EXO Treatment

5 μ L of Treated PCR Reaction



II Multiplex ASPE

5 μ L of ASPE Reaction



III Bead Hybridization

Filtration and Wash

**IV Data Acquisition
on xMAP**



Data Analysis

1 hr


15 min

45 min



5 min





Tag-It Multiplex (16-plex) PCR




CFTR 40+4: Net Signals (MFI)




CFTR 40+4: Allelic Ratios




Tag-It Specificity



Correct Calls as a Function of S:N



Tag-It Data Acquisition





TDAS Summary View

TDAS (version 4.01) - Summary view

File View Sample Variation Help

File name: C:\My Sessions\CFTR ExampleData1\Output.csv Date and time created: 06/25/2003,10:50:40 AM
Session name: ExampleData1 Created by: CF

Lo...	Sample	Mu. alleles det...	Mu. and wt all...	G85E	394deITT	R117H	Y122X	I148
A1	Sample 1			WT	WT	WT	WT	V
B1	Sample 2			WT	WT	WT	WT	V
C1	Sample 3		394deITT, dF508	WT	HET	WT	WT	V
D1	Sample 4	R347H	R117H	WT	WT	HET	WT	V
E1	Sample 5			WT	WT	WT	WT	V
F1	Sample 6	dI507, F508C		WT	WT	WT	WT	V
G1	Sample 7	2183AA>G, 21...		WT	WT	WT	WT	V
H1	Sample 8		Y1092X-C>A	WT	WT	WT	WT	V
A2	Sample 9			WT	WT	WT	WT	V

Analyzed using: **CFTR**

- Only wild-type allele detected (or Call Hidden)
- Mutant allele(s) detected
- Wild-type and mutant alleles detected
- Sample has at least one mutant allele detected
- Call(s) not made
- Signal significant
- Control sample

This data file has been analyzed using the Tm Bioscience CFTR analysis module. The possible genetic calls are WT (only wild-type allele detected), Mu D (mutant allele detected), HET (heterozygous, both wild-type and mutant alleles detected or two different mutant alleles detected), ND (no mutant allele detected for the I506V, I507V, F508C variation), CH (call hidden), or No Call (no call possible). For more information about the genetic calls, the messages, and the software in general, please read the Help or the manual.
CAUTION: As with any hybridization-based assay, underlying polymorphisms in primer-binding regions can affect the alleles being detected

For Help, press F1

NUM

Genetic calls are displayed for each variation, for all samples.



Detailed Sample View

Complete data for sample D1 (Sample 4) on session ExampleData1

Variation	Call	Raw Signal (MFI)		Background (MFI)		Net Signal (MFI)		Allelic Ratio		AR Thresholds		
		Wt Allele	Mu Allele	Wt Allele	Mu Allele	Wt Allele	Mu Allele	Wt Allele	Mu Allele	WT Call	Wt Present	Mu Present
G85E	WT	2224.5	241.0	36.5	33.5	2188.0	207.5	0.91	0.09	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
394delTT	WT	1380.5	63.5	35.0	47.5	1345.5	16.0	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
R117H	HET	601.0	757.0	23.0	30.5	578.0	726.5	0.44	0.56	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
Y122X	WT	1028.0	54.0	37.0	45.0	991.0	9.0	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
I148T	WT	1430.0	69.0	25.0	37.0	1405.0	32.0	0.98	0.02	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
621+1G>T	WT	1548.5	52.0	37.0	49.0	1511.5	3.0	1.00	0.00	0.85 - 1.00	0.25 - 1.00	0.30 - 1.00
711+1G>T	WT	1492.0	42.0	27.0	30.0	1465.0	12.0	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
1078delT	WT	1701.0	93.5	29.0	36.0	1672.0	57.5	0.97	0.03	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
R334W	WT	955.0	75.0	43.0	30.5	912.0	44.5	0.95	0.05	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
R347P	Mu D (H)	56.0	67.0	33.0	40.0	23.0	27.0	0.01	0.02	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
└ R347H			1673.5		32.0		1641.5		0.97			0.30 - 1.00
A455E	WT	1924.0	49.0	29.0	39.0	1895.0	10.0	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00

Analyzed using: **CFTR**

- Only wild-type allele detected (or Call Hidden)
- Mutant allele(s) detected
- Wild-type and mutant alleles detected
- Sample has at least one mutant allele detected
- Call(s) not made
- Signal significant
- Control sample

This data file has been analyzed using the Tm Bioscience CFTR analysis module. The possible genetic calls are WT (only wild-type allele detected), Mu D (mutant allele detected), HET (heterozygous, both wild-type and mutant alleles detected or two different mutant alleles detected), ND (no mutant allele detected for the I506V, I507V, F508C variation), CH (call hidden), or No Call (no call possible). For more information about the genetic calls, the messages, and the software in general, please read the Help or the manual.

Print Help Close Close all sample data windows

All data and genetic calls are displayed for the selected sample.



Detailed Variation View

Complete data for variation R117H on session ExampleData1

Sample		Call	Raw Signal (MFI)		Background (MFI)		Net Signal (MFI)		Allelic Ratio		AR Thresholds		
Location	Name		Wt Allele	Mu Allele	Wt Allele	Mu Allele	Wt Allele	Mu Allele	Wt Allele	Mu Allele	WT Call	Wt Present	Mu Present
A1	Sample 1	WT	1037.5	57.0	23.0	30.5	1014.5	26.5	0.97	0.03	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
B1	Sample 2	WT	3000.0	65.0	23.0	30.5	2977.0	34.5	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
C1	Sample 3	WT	1757.0	60.0	23.0	30.5	1734.0	29.5	0.98	0.02	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
D1	Sample 4	HET	601.0	757.0	23.0	30.5	578.0	726.5	0.44	0.56	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
E1	Sample 5	WT	3162.0	64.5	23.0	30.5	3139.0	34.0	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
F1	Sample 6	WT	1757.0	60.0	23.0	30.5	1734.0	29.5	0.98	0.02	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
G1	Sample 7	WT	1757.0	60.0	23.0	30.5	1734.0	29.5	0.98	0.02	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
H1	Sample 8	WT	1757.0	60.0	23.0	30.5	1734.0	29.5	0.98	0.02	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
A2	Sample 9	WT	2613.5	49.0	23.0	30.5	2590.5	18.5	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
B2	Sample 10	WT	2559.0	50.0	23.0	30.5	2536.0	19.5	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
C2	Sample 11	WT	2948.0	58.5	23.0	30.5	2925.0	28.0	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
D2	Sample 12	WT	3007.5	59.0	23.0	30.5	2984.5	28.5	0.99	0.01	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
E2	Sample 13	HET	1253.0	1531.5	23.0	30.5	1230.0	1501.0	0.45	0.55	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00
F2	Sample 14	Mu D	32.0	2400.0	23.0	30.5	9.0	2369.5	0.00	1.00	0.85 - 1.00	0.30 - 1.00	0.30 - 1.00

Analyzed using: **CFTR**

- Only wild-type allele detected (or Call Hidden)
- Mutant allele(s) detected
- Wild-type and mutant alleles detected
- Sample has at least one mutant allele detected
- Call(s) not made
- Signal significant
- Control sample

This data file has been analyzed using the Tm Bioscience CFTR analysis module. The possible genetic calls are WT (only wild-type allele detected), Mu D (mutant allele detected), HET (heterozygous, both wild-type and mutant alleles detected or two different mutant alleles detected), ND (no mutant allele detected for the I506V, I507V, F508C variation), CH (call hidden), or No Call (no call possible). For more information about the genetic calls, the messages, and the software in general, please read the Help or the manual.

Print Help Close Close all variation data windows

All data and genetic calls are displayed for the selected variation.



Ashkenazi Jewish Panel

- Bloom Syndrome
- Canavan Disease
- Familial Dysautonomia
- Fanconi Anemia Group C
- Gaucher Disease
- Mucopolysaccharidosis Type IV (ML-IV)
- Niemann-Pick Disease
- Tay-Sachs disease

Additional Products

- PGx Panels (2D6, 2C9, 2C19)
- Fragile X
- Gene Dosage (aneuploidy)
- ID (MRSA, HPV)

Tag-It™ Advantage

- Multiplexed—entire assay (86 alleles) tested in one well
- Complete genotyping—analyzes both wildtype & mutant at each locus in each assay
- Accuracy—S:N of 20:1 to 50:1 for each allele
- Automation—Up to 1200 samples/day/technician in a 96-well format
- Menu—CF40, Jewish panel (8 diseases), P450 testing (2D6, 2C9, 2C19), others



Tag-It™ Advantage

- No more reflexing—reflex tests included in each assay, but results only shown if needed
- Full data analysis software provided; makes summary calls; can easily interface with LIS
- Low-cost instrument with ability to process 96 samples/40 minutes





Putting the human genome to work™

www.tmbioscience.com