

UKNEQAS - Molecular Diagnosis of Haematological Malignancies

BCR-ABL Quantitation Programme

Jane Holden



MDHM Programme

- ◆ 110 participants registered
- ◆ Organised into 5 sub-programmes:
 - ◆ IgH/TCR clonality testing
 - ◆ *JAK2* V617F status
 - ◆ BCR-ABL and AML translocation identification (including t(8;21), t(15;17) and inv(16))
 - ◆ BCR-ABL quantitation
 - ◆ Post SCT chimerism monitoring
- ◆ 2 trials issued per year per sub-programme



BCR-ABL Quantitation Programme

- ◆ 63 participants
- ◆ K562 (M-BCR) cell-line distributed for BCR-ABL quantitative analysis – lyophilised to ensure RNA stability
- ◆ 2 trials issued




BCR-ABLQ Trial 1


- ◆ Issued to 39 participants; only 24 returned results
- ◆ 1 sample issued – 1% dilution of K562 cells in HEL cells (lyophilised)
- ◆ Participants asked to perform quantitative analysis and report M-BCR level



BCR-ABLQ Trial 1



BCR-ABLQ Trial 1




BCR-ABLQ Trial 2

- ◆ Issued to 63 participants; 55 returned results
- ◆ 3 samples issued (lyophilised)
 - 100% K562 sample (diagnostic sample)
 - ~0.4% dilution of K562 in HEL (follow-up sample)
 - 100% HEL sample (negative sample)
- ◆ Participants asked to perform quantitative analysis and report Log reduction in M-BCR at follow-up



BCR-ABLQ Trial 2

Graph to show the distribution of Log reduction results for MDHM BCR-ABLQ



BCR-ABL Quantitation Results

- ◆ Variation in quantitative level due to:
 - Different in-house protocols/kits
 - Different control genes
 - Different material for standard dilutions



Method Parameter	Number of Participants
Control Gene:	
ABL	42
G6PDH	6
Beta-2 microglobulin	3
β-GUS	3
BCR	1
Material for Standard Dilutions:	
Commercial Plasmids (Ipsogen)	32
In-house Plasmids	13
Serial dilution from K562 DNA	3
No standards used	4
No info provided	3


BCR-ABL Quantitation Results

- ◆ Variation in quantitative level due to:
 - Different in-house protocols/kits
 - Different control genes
 - Different material for standard dilutions
 - Different instruments and analysis software



Instrument	Number of Participants
ABI 5700	1
ABI 7000	6
ABI 7300	2
ABI 7500	10
ABI 7700	4
ABI 7900	8
Roche Light Cycler	15
Rotor-Gene	4
Cepheid Smartcycler	1
BioRAD Icyler	1
Stratagene MX3000p	2

BCR-ABLQ Trial 2




BCR-ABL Quantitation Results

- ◆ Variation in quantitative level due to:
 - Different in-house protocols/kits
 - Different control genes
 - Different material for standard dilutions
 - Different instruments and analysis software
 - RNA/cDNA quality



BCR-ABLQ Trial 2



BCR-ABL Quantitation Results

- ◆ Control gene Ct values vary greatly indicating variation in the quality of RNA/cDNA:
 - ABL 19.93-32.08
 - G6PDH 19.12-35.95
 - β -GUS 22.26-25.76
 - B2M 24.69-25.53
 - BCR 34.98
- ◆ Several labs perform RQ-PCR with only one or two replicates
- ◆ Quantitative M-BCR level is not expressed in a standard way



BCR-ABL Q Future Work

- ◆ Decrease M-BCR level in trial samples
- ◆ Introduce scoring system – enables consistently poor performers to be identified
- ◆ Introduction of standardisation and guidelines



Acknowledgements

UKNEQAS LI

Viv Granger

Liam Whitby

Alison Whitby

Joanne Antcliffe

Sarah Sardinha

Dr David Barnett

Prof John Reilly

Sheffield Molecular Genetics
Service

University of Sheffield Medical
School (Haematology dept)

UKNEQAS Genetics SAG

jlholden@btconnect.com

