

**Economical liquid phase dose assays in microplates:
proofs of principle for sex testing and haptoglobin
duplicon genotyping**

Philip Guthrie, Bristol Genetic Epidemiology Lab.

**>1400 variable loci at last count
(Sanger Centre)**


Segmental duplication content of the human genome

Chromosome	% Intra	% Inter
1	2.1	1.2
2	2.0	1.4
3	1.1	0.8
4	1.4	1.1
5	1.9	0.8
6	1.4	0.4
7	5.5	1.7
8	1.6	0.8
9	5.5	3.6
10	3.9	1.0
11	2.6	1.2
12	1.4	0.7
13	0.8	1.1
14	0.5	0.8
15	4.6	2.4
16	6.0	2.7
17	6.0	0.8
18	0.7	0.9
19	4.5	1.2
20	0.9	1.4
21	1.1	2.9
22	3.6	2.9
X	1.8	2.0
Y	10.2	6.0

Defensins

- **Important part of immune system**
- **Small peptides**
- **Hexamers – ion channels in pathogen membrane**
- **Permeabilisation → cell death**
- **DEFA1: 3 - 9 copies (majority 5)**
- **DEFA3: 0 – 4 copies (majority 1 or 2)**


CNV in Defensins



Cytochrome P450 2D6 (CYP2D6)

- Involved in metabolism of up to 20% of all drugs in clinical use
- 0 copies – ‘Poor’ metabolisers
- 1 copy – ‘Intermediate’
- 2 copies – ‘Extensive’
- 3 copies – ‘Ultra-rapid’
- Low copy number – high drug plasma concentrations
- High copy number – rapid drug degradation & therapeutic failure

Haptoglobin



2/2 genotype:

- reduced vitamin C levels
- disease associations e.g. diabetic vascular complications

1/1



1/2




2/2






- 4 primers – 4 different properties
- – original template ratios??
- Numerous possible confounding factors –
 - DNA quantity, DNA quality
 - Different fragment sizes
 - Protein
 - EDTA
 - Mg
 - pH
 - log effect of different amp. efficiencies over e.g. 35 cycles



- Better
 - 1 forward universal primer
 - 1 reverse universal primer


MADGE gel assay for Haptoglobin junction null allele



Roche LightTyper




Differential melting of A/T- or G/C-rich fragments




Target gene:- ~49% G/C

Reference gene:- ~58% G/C


Trial of Target/Reference fragment differential melting




X chromosome complement, by melting profile




Haptoglobin duplication - null allele homozygous




Haptoglobin duplication allele - heterozygous




Haptoglobin duplication allele - homozygous




Haptoglobin duplication – composite profiles




Hp Junction/Reference fragment – area under curve ratios




384-well microplate fluorescing in LightTyper




BWHHS 384-well microplate – Hp Jcn/Ref peak height ratios




BWHHS 384-well microplate – Hp Jcn/Ref ratios



Hp Jcn/Ref ratios closeup




Hp Jcn/Ref ratio gradients






- 4 primers – 4 different properties
- – original template ratios??
- Numerous possible confounding factors –
 - DNA quantity, DNA quality
 - Different fragment sizes
 - Protein
 - EDTA
 - Mg
 - pH
 - log effect of different amp. efficiencies over e.g. 35 cycles



- Better
 - 1 forward universal primer
 - 1 reverse universal primer



- Best
 - developments currently being patented
 - sorry, can't disclose yet

Acknowledgments

Prof. Ian Day

Dr. Reza Abdollahi

Dr. Tom Gaunt

Dr. Xiao-he Chen

Prof. Shah Ebrahim

Dr. Debbie Lawlor

Prof. George Davey-Smith

Ms. Rita Patel

BGEL startup was supported by the University of Bristol.

BWHHS is funded by the UK Department of Health.

Genetic studies in BWHHS have been supported by BHF.

Work in the laboratory of INMD, TRG and XC is supported by UK MRC and BHF.